- (GB)** 0 800 731 1792
- (IE)** 1 800 509 448
- (FR)** 0 800 944 803  
(appel gratuit depuis un poste fixe)
- (BE)** 0 800 14 592
- (LU)** 4 00 50 51
- (TR)** 0 800 261 63 65
- (RU)** 8-800-200-20-20
- (UA)** +38 044 428 65 05
- (MA)** 0800 100 70 75

**Internet:**

[www.oralb.com](http://www.oralb.com)

[www.braun.com](http://www.braun.com)

[www.service.braun.com](http://www.service.braun.com)

Braun GmbH

Frankfurter Straße 145

61476 Kronberg/Germany

**Type 3757**

D12.513

D12.513W

D12.513S

99325376/VII-11

GB/F/TR/RO/MD/RUS/UA/KZ/Arab

2

## English

Welcome to Oral-B!

Please read the use instructions carefully before first use.

### IMPORTANT

- Periodically check the cord for damage. If the cord is damaged, take the charging unit to an Oral-B Braun Service Centre. A damaged or non-functioning unit must no longer be used.
- This product is not intended for use by children under age three.
- This appliance is not intended for use by children or persons with reduced physical, sensory or mental capabilities, unless they are supervised by a person responsible for their safety. In general, we recommend that you keep the appliance out of reach of children.
- Children should be supervised to ensure that they do not play with the appliance.
- If the product is dropped, the brush head should be replaced before the next use even if no damage is visible.
- Do not place or store the charger where it can fall or be pulled into a tub or sink.
- Do not place the charger in water or other liquid.
- Do not reach for a charger that has fallen into water. Unplug immediately.
- Do not modify or repair the product. This may cause fire, electric shock or injury. Consult your dealer for repairs or contact an Oral-B Service Centre. Do not disassemble the product except when disposing of the battery.
- When taking out the battery for disposal, use caution not to short the positive (+) and negative (-) terminals.
- Do not insert any object into any opening of the appliance/charger.
- Do not touch the power plug with wet hands. This can cause an electric shock.
- When unplugging, always hold the power plug instead of the cord.
- Use this product only for its intended use as described in this manual. Do not use attachments which are not recommended by the manufacturer.
- If you are undergoing treatment for any oral care condition, consult your dental professional prior to use.

### Description

- A Brush head (depending on model/ interchangeable)
- B On/off switch
- C Handle
- D Charging unit

### Connecting and charging

Your toothbrush has a waterproof handle, is electrically safe and designed for use in the bathroom.

- Plug the charging unit (D) into an electrical outlet and place the handle (C) on there. For voltage specifications refer to the bottom of the charging unit. A full charge takes 16 hours and enables approximately one week of regular brushing (twice a day, 2 minutes) (1).
- For everyday use, the handle can be stored on the plugged-in charging unit to maintain the toothbrush at full power; there is no risk of over-charging. However for environmental reasons, Oral-B recommends to unplug the charging unit until the next recharge is required.

### Battery maintenance

To maintain the maximum capacity of the rechargeable battery, unplug the charging unit and discharge the handle by regular use at least every 6 months.

### Using the toothbrush

Wet brush head and apply toothpaste; you can use any kind of toothpaste. To avoid splashing, guide the brush head to your teeth before switching the appliance on (B) (2).

Start brushing the outsides, and then the insides and finally the chewing surfaces (3). Turn off the handle by pressing the on/off switch.

#### Timer (depending on model)

An incorporated timer memorises the elapsed brushing time. A short stuttering sound indicates the end of the professionally recommended 2-minute brushing time.

### Individual brush heads and accessories

Oral-B offers you a variety of brush heads which are available in stores that sell Oral-B products:

#### The Oral-B Precision Clean brush head

cups each tooth with a curved bristle trim and reaches deep between teeth at the same time.

#### The Oral-B Sensitive brush head

an extra soft bristle texture that is gentle on teeth and gums.

#### Oral-B 3DWhite brush head

features a dentist-inspired polishing cup to remove surface stains from coffee, tea and tobacco for effective brightening. We do not recommend using the Oral-B 3DWhite brush head with braces.

Dentists recommend replacing any toothbrush every 3 months. Therefore brush heads feature light blue INDICATOR® bristles to help you monitor brush head replacement need. With thorough brushing, twice a day for two minutes, the blue color will fade halfway down approximately in a 3-month period. If the bristles splay before the colour recedes, you are exerting too much pressure on teeth and gums (3a–c).

During the first days of using the toothbrush, your gums may bleed slightly. In general, bleeding should stop after a few days. Should it persist after two weeks, please consult your dentist or dental hygienist.

### Cleaning

After use, rinse the brush head for several seconds under running water with the handle switched on. Then switch off the handle and take off the brush head. Rinse both parts separately under running water and wipe them dry. From time to time, also clean the charging unit with a damp cloth (4).

Subject to change without notice.

This product contains rechargeable batteries. In the interest of protecting the environment, please do not dispose of the product in the household waste at the end of its useful life. Disposal can take place at an Oral-B Braun Service Centre or at appropriate collection points provided in your country.

## Garantie

We grant 2 years guarantee on the product commencing on the date of purchase. Within the guarantee period we will eliminate, free of charge, any defects in the appliance resulting from faults in materials or workmanship, either by repairing or replacing the complete appliance as we may choose. This guarantee extends to every country where this appliance is supplied by Braun or its appointed distributor.

This guarantee does not cover: damage due to improper use, normal wear or use, especially regarding the brush heads, as well as defects that have a negligible effect on the value or operation of the appliance. The guarantee becomes void if repairs are undertaken by unauthorized persons and if original Braun parts are not used.

To obtain service within the guarantee period, hand in or send the complete appliance with your sales receipt to an authorized Oral-B Braun Customer Service Centre.

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